

OR1E2-olfactory receptor, family 1, subfamily E, member 2 Gene

Reference: PTXBC121105

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OR17-135; OR17-136; OR17-93; OR1E4; OR1E7; OST529; olfactory receptor 1E2; olfactory receptor 17-93/17-135/17-136; olfactory receptor 1E4; olfactory receptor OR17-17; olfactory receptor OR17-3; olfactory receptor, family 1, subfamily E, member 4; olfactor		

For research use only.

GENE DESCRIPTION	olfactory receptor, family 1, subfamily E, member 2
GENE ID	8388
ORIGIN SPECIES	Human
ACCESSIONS	BC121105
VECTOR	pDONR223
NCBI SYMBOL	OR1E2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgggacaaaatcaaccagcatctcagacttctgtctcctggcctgccatccaaccagagcagcaaaacctgcgtatgccctgttcttgccatgtatctaccacctctggggaacctctcatcattgtctcattcgactggactcccatctccacacgectgtgtattgtttct
cagcaactgtctctctgacctctgttttctcagtcacaatgcccaatgctgcagaacatgcagaaccaagaccatccatccctatgcagactgctgacccaatgtacttcttctgtatttttcggatctagagagcttctccttggccatggcctatgaccgctatgtggccatc
tgcctcccatgcaactacaccgcatctgcttctcctgactacaccgcatcatgagcccatgctctctcctggtggcgctgtctcctgggtgctgaccaccttccatgcatgttacacactttactcatggccagggtgtgttttgcagacaatgtgatccccactttttctgtgat
gtctgtctgtgaagctggcctgtctgacactcagtttaagtgaatgggtgatattatcatggggggtcattctgtcatccattctactcatccttgggtctatgcaagaattgtctcctcatcctcaaggctccttcttaagggtatctgaaggccttctactgtgtgctccacct
ctctgtgtgtcactgttctatggaccgtattgtctactatgcccacagetaatagttctactctaaaggacactgtcatgctgatgtacactgtgtgaccccatgctgacccctcatctacagcctgaggacagagacatgaaggagccctggaagggtcatttga
aggaaaaatcccttctctatga