

LOC340017-hypothetical protein LOC340017 Gene

Reference: PTXBC121112

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC340017		
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC340017
GENE ID	340017
ORIGIN SPECIES	Human
ACCESSIONS	BC121112
VECTOR	pDONR223
NCBI SYMBOL	LOC340017
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggttctctggtgtctcaagcttcggatgccactgattctcagtgccagccgtggaagctgctgtgtacatttggccagccacagccttcagggagctgggtcccatgctggtgcctatagctccccgcccccaccacagccagcttgccatggccatgcacagtggccagacctcat
gctcactcatgacaccttcgccactccgacctggcccgcccttggcaggtatgggattcagggatccagccccagccgattgcagcctgccagggcagtggtggaatgagacctgtggcccaagcaaaactcaggcaaaaggcaccactggccacatag