



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

UBAC2-UBA domain containing 2 Gene

Reference: PTXBC121139

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

PHGDHL1; ubiquitin-associated domain-containing protein 2; RP11-178C10.1; UBA domain-containing protein 2; phosphoglycerate dehydrogenase like 1; phosphoglycerate dehydrogenase-like protein 1; UBA domain containing 2

For research use only.

GENE DESCRIPTION	UBA domain containing 2
GENE ID	337867
ORIGIN SPECIES	Human
ACCESSIONS	BC121139
VECTOR	pDONR223
NCBI SYMBOL	UBAC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagatattcttggcatcactgcagctagtaatttgcctctggattctcggcactctgtttgctctgtttgaccatttactgccataccagagctccaagtggcacaattctgggtccgttgcctatcacaacaagacattgatttataattgggactgcagctttccacctctggttcctatcatctggattgtagccataaaggaggacttattcgcgctgtctgtctacgacagcaaaatgtccagggtgcacaggtctctgcatcccgctgatggcaaaatttcttctggacactgaaccacattctcttcttcagaaacccacagcgaagccagaattgggatggagaccagctgtgga catccagacagacagagaagaagcagctgctgtggccgcgacgtctgattgtctcagttgttcacaaagggtggcgcgacagacagcagcagcgaaggaatgactcaattggaaatctgtttctttacgtcagcgcacaaacgtaaactcagggcggctgcggagctctgagccagcagc gccccctagagaatttctgtagaacacgctgtgcccgccctatgagatggattgttccagaaaggatggcgcgcctgagacgtctcaaaacgaactcaatgtgccacattctctctctgacagcactga