

LOC285045-hypothetical LOC285045 Gene

Reference: PTXBC121144

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC285045
GENE ID	285045
ORIGIN SPECIES	Human
ACCESSIONS	BC121144
VECTOR	pDONR223
NCBI SYMBOL	LOC285045
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagaaaacagtggcttcaagggttaacaactgtcttcgattgcacagctaacaagagctcactgtccaaaggaggtctctctccccgtccacatgctgaggcttgcctccagagtgaatgtgtgactggaagaaccacacagggccctgagagctgtcagcctgccagcttctgtgagatccctgggtccgagtagccaagctggaggagcaccatttgcctactgactcttgggcaacgtgcagcctgcatccagccggcagaaggcaggtggtatcacagctgagtttacagccacggaaggaaaccactgttagcctttag