

LOC284276-hypothetical LOC284276 Gene

Reference: PTXBC121150

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC284276
GENE ID	284276
ORIGIN SPECIES	Human
ACCESSIONS	BC121150
VECTOR	pDONR223
NCBI SYMBOL	LOC284276
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaccgggttggtctgtttccctggaagcactgctgggaaccaagccagccctgcagggtctctctgctgcccagccagcagcctctcttgcgctgcagcccgctgctgcaccgggttcatcggctgaagtcagctgctctcttggcctctgcacctca
gcacaggctactgctctgagggccgcccagccgctgctctctcacaacctccatgatcgtggccctgaggtcaccttccattggcgatgaacaatcccatgtcactcccatgggtgtgccacgggaagctgctcatgggcagggagggaactgttcaagctggttttccagcc
cctaccacagcgtgggtgtgtaggtggacaagaaagtcattacgtatgtaaatga