



Reference: PTXBC121178

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AI1H; integrin beta-6; integrin, beta 6; integrin subunit beta 6		
For research use only.			

GENE DESCRIPTION	integrin, beta 6
GENE ID	3694
ORIGIN SPECIES	Human
ACCESSIONS	BC121178
VECTOR	pDONR223
NCBI SYMBOL	ITGB6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr