



Reference: PTXBC121179

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DMRV; GLCNE; IBM2; Uae1; bifunctional UDP-N-acetylglucosamine 2-epimerase/N-acetylmannosamine kinase; N-acylmannosamine kinase; UDP-GlcNAc-2-epimerase/ManAc kinase; UDP-N-acetylglucosamine 2-epimerase/N-acetylmannosamine kinase; glucosamine (UDP-N-acetyl)		

For research use only.

GENE DESCRIPTION	glucosamine (UDP-N-acetyl)-2-epimerase/N-acetylmannosamine kinase
GENE ID	10020
ORIGIN SPECIES	Human
ACCESSIONS	BC121179
VECTOR	pDONR223
NCBI SYMBOL	GNE
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]