

TAF13-TAF13 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 18kDa Gene

Reference: PTXBC121180

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TAF13 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 18kDa; TAF(II)18; TAF2K; TAFII-18; TAFII18; transcription initiation factor TFIID subunit 13; TATA box binding protein (TBP)-associated factor, RNA polymerase II, K, 18kD; transcri		

For research use only.

GENE DESCRIPTION	TAF13 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 18kDa
GENE ID	6884
ORIGIN SPECIES	Human
ACCESSIONS	BC121180
VECTOR	pDONR223
NCBI SYMBOL	TAF13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagatgaggagaagacccacgtttgaggagaagaatgaagaattggaggaggtgcagaaggtgcacagggtaaaagaagagacgttttctaaagaattgcgatgatgatgtatggcttggggatgaccagaatccttatactgagtcagtgatattctgaagatcttgca
tagagtttatcactgaaatgactcacaaggcaatgtcaattggaagacaaggcagtagacaagttgaagatactgtcttctgattcgaaggacccaaggaagttgccagggttaagactgcttactatgaatgaagaattgaacgagctagaaaagcatttgatgaagcaaatatggat
cttga