



Reference: PTXBC121800

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BT3.1; BTF5; BTN3.1; CD277; butyrophilin subfamily 3 member A1; dJ45P21.3 (butyrophilin, subfamily 3, member A1)		

For research use only.

GENE DESCRIPTION	butyrophilin, subfamily 3, member A1
GENE ID	11119
ORIGIN SPECIES	Human
ACCESSIONS	BC121800
VECTOR	pDONR223
NCBI SYMBOL	BTN3A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]