

AMTN-amelotin Gene

Reference: PTXBC121817

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UNQ689; amelotin; RSTI689		
For research use only.			

GENE DESCRIPTION	amelotin
GENE ID	401138
ORIGIN SPECIES	Human
ACCESSIONS	BC121817
VECTOR	pDONR223
NCBI SYMBOL	AMTN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagggtacgattctactgtttgtctcttaggatcaactcgtcattaccacagctcaaacctgttgggactccctccacaaaactggtccggatcaggggaacactacaaaccaacagcagtcgaatcaggctcttctctttaagtctgataccattaacacagatgctcacactgg
ggccagatctgcatctgttaaatcctgctgcagggaatgacacctgtgtaccagaccacccattgacctgggaggttgaatgtacaacagcaactgcacccacatgtgtaccaattttgtcacacaactggagccaggggactatcctaagctcagagggaattgccacaaatcttcac
gagcctcatcatccattctgttcccgaggagcatctgcccaccagtcaggcaggggctaatccagatgtccaggatggaagccttcagcagaggagcaggtgtaaatctgccacccaggggaacccagcagggcgctcccaactccagtgccacagatgacgacttgc
gtgaccacccctgcaggcatccaaaggagcacacatgccattgagggaagccaccagaaatcagaatggaattcagtaa