

FLJ42957-FLJ42957 protein Gene

Reference: PTXBC121822

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NCRNA00173; long intergenic non-protein coding RNA 173		
For research use only.			

GENE DESCRIPTION	FLJ42957 protein
GENE ID	400077
ORIGIN SPECIES	Human
ACCESSIONS	BC121822
VECTOR	pDONR223
NCBI SYMBOL	FLJ42957
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcttctcccggtacagcacgatgattacagctgtgtttgttcaacagtcgtgtacaactgacagtgccatcttactgcctggctcaggtcaggttacttaaggcttattatgtgttacgaagggcagcacaggaaaaggacaaggggtgtgtcagggtggcactgtgttaaaaa
gtggcggtgcaagggccgcatccggcgagccgtgcaacctcagccctggcccttacctcgcagcctctccagcatccagctaccagactccaagggccagcgagagccagctctcgggtacgtggagctccacaggtccagaatcggtggatcagagttcaaatc
tggttctgctactgtctaattgcgtgctgcagggactcaatctctcatctgggaatgggagtaa