

C8orf77-chromosome 8 open reading frame 77 Gene

Reference: PTXBC122544

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C8orf77; ZNF252P antisense RNA 1		
For research use only.			

GENE DESCRIPTION	chromosome 8 open reading frame 77
GENE ID	286103
ORIGIN SPECIES	Human
ACCESSIONS	BC122544
VECTOR	pDONR223
NCBI SYMBOL	C8orf77
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccggcagtcctgtctctccgggttacgtccctccagcccttggcggagctcgcggcgagaggcgaggtgcagaggtgcggcgccggcggtgcggctgcgaaggaggaccccgacaccgagcgttctgcggacgtcctcaaccgggggtgttctcgtgctctg
gacccggccctctccacctaggacgtctcgcgggtgcacgtgccccctctcgcggcgactttcccccaacctccaaagtatgatcaagctcagaatccagctcgcggccgctgtgcacctagcacttcccacctgtagtctttgacctcatgccccctctagaaggctgtccatcgtta
aatgagggaagccgagggtgcagagaggtttcagacctattgcgggtggagctggagttcgagaaccaaccctcggctcgtgagacgcgcttgcggaacggcgccggcgccgggggtgaagcggagtgaggggagagggccaggtgcgacccgggcaggtgaggagcacggggccg
gaggggaggacttacgagaatggaaaggaaagcagccgtctcagttgggacagtgaggagcattaatgtcagagaatgagcagttgtgggaagagccatag