



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ZNF260-zinc finger protein 260 Gene

Reference: PTXBC122549

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

OZRF1; PEX1; ZFP260; zinc finger protein 260; zfp-260

For research use only.

GENE DESCRIPTION	zinc finger protein 260
GENE ID	339324
ORIGIN SPECIES	Human
ACCESSIONS	BC122549
VECTOR	pDONR223
NCBI SYMBOL	ZNF260
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]