

LOC284749-hypothetical protein LOC284749 Gene

Reference: PTXBC122555

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	long intergenic non-protein coding RNA 494		
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC284749
GENE ID	284749
ORIGIN SPECIES	Human
ACCESSIONS	BC122555
VECTOR	pDONR223
NCBI SYMBOL	LOC284749
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagcatggaaatgctgtgtgagtgagctaagcagccgttctcaggagaaaacagacaaatcctcctgcctcctgcaaacagccttggcctgttttatgtcatctttgtgagctccacgggtcacccctcctgttggtttgtgatactgacatttcaggcactggcacagatacatt
ggccttctggtggtgatttcacggacggagtgagcagtgctggttggaggagccgtggcgccacagggcctgtgctcagctgatcaacgccttgatagccaggaaaggcctgtcctaccccaaggccctcgagcactgtttgaaattgctcttgcctcttggaaatggattcct
cagtacatttcagactcgtgtgttga