

FLJ31222-hypothetical LOC388387 Gene

Reference: PTXBC122868

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC388387
GENE ID	388387
ORIGIN SPECIES	Human
ACCESSIONS	BC122868
VECTOR	pDONR223
NCBI SYMBOL	FLJ31222
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaactggccagatcaactgctgtcaccttctaaagagagggttttgcctgttgcaggctggctcgaactcctggactcaagtgtatgctgctcagcctccaaagtgtggaattataggaccagacggggaggagggaatgccagagaagctgcagactctccggccac
ggagcctgagtgagcagaggtggccactgtcagccacaggtcaagtgtcagggtgggacggggcaggctccaccagaaggcggggcaggctccaccagaaggcggggccagaccacagggcaagggtggagatgagctgtgtgcacccccaggaagttaccagattgc
cccagagcttctgacctttggatctctctgggtgcatctctcagggtcccaatgtcttcttgaagggtggaaggaaatgtcccatcagagcctccaacaggccacctgggggtctctctga