

ZNF586-zinc finger protein 586 Gene

Reference: PTXBC124558

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein 586		
For research use only.			

GENE DESCRIPTION	zinc finger protein 586
GENE ID	54807
ORIGIN SPECIES	Human
ACCESSIONS	BC124558
VECTOR	pDONR223
NCBI SYMBOL	ZNF586
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcagcagccgctctgagggcgctcaggtgtgtggcatggagggaagatgagcagcagccttctaagcagagcagctgtatatacatatatacaagaccaggagggtcatagtgagaaaggccttatgagtggtgggaatataggaaattttaagaacaagtcctgcc
tcactgaaccgaagaagatcacaaacacagggaatgtcgactggagaaaggccttatgagtcagcaaatatgggaattattcaccaaaagcctacactccatattcatgagagatttcatactgggcaaaagacctatgagtcagtgagtggtgaaatcattcaccaaagctctc
actcttgacgctgacacacttcacactagagaaaggccttatgagtgattgaatgtgggaaggcctttgctgaaaagtcagctcattaccacaggaaggtcactctggagcaaaagcgttatgaatgcaatgaatgtgggaagtcctttgctatatacatctagtcattaaacacaggag
gattcacactggagaggccttatgagtcagtgaaatgtgggagatcctttgctgaaaactccagtccttattaaacactga