

SGCZ-sarcoglycan zeta Gene

Reference: PTXBC125037

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZSG1; zeta-sarcoglycan; zeta-SG; sarcoglycan zeta		
For research use only.			

GENE DESCRIPTION	sarcoglycan zeta
GENE ID	137868
ORIGIN SPECIES	Human
ACCESSIONS	BC125037
VECTOR	pDONR223
NCBI SYMBOL	SGCZ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacacgagacaatacactagcaacccaacagaataaactgccaaggactgagaatgcacaactttaccagtgagggaatttatggatggcgaagagggtcttatactctttgtccttctgtgttggtaccatgatagtttaactagccatgacaatatggatattgaaagtatgaattt
cactgttgatagtcgctgtgtttacagctgtgacaggaatgtcacagtgaatgcaagaaatcacatggggcagtttaaccggacagctgaccataggagctgatgtgtggaagctcagtgtaaaagattgaaagtgagagccagtgaaatggcaggggtctgttttctgcagatgaagatg
agattaccattggggctgaaaagctgaaagtacaggcactgaaggagccgtatttggcactctgtggagacgcccacatcagagcagagagccatcccaagatctcaggcttgatcaccaccagatccttgatcatggaagctccccgtgggtccagtgagtgctgctgcaggag
actcaaggccactgcaggaaggagctccatctgcaatctacagaaggaggagatatttttaaatgcagagacaatcaagctgggaatctaccaactggctcctctctcttcttccaccagctctcaagttctcgcagacagtgatgaactctgctctgccccaatggcaactttacc
tttctccagcaggagtaggttccactgtcagtcagtagcaacatctgctgtggagctga