

BAG2-BCL2-associated athanogene 2 Gene

Reference: PTXBC125039

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BAG-2; dJ417I1.2; BAG family molecular chaperone regulator 2; bcl-2-associated athanogene 2; dJ417I1.2 (BAG-family molecular chaperone regulator 2); BCL2 associated athanogene 2		

For research use only.

GENE DESCRIPTION	BCL2-associated athanogene 2
GENE ID	9532
ORIGIN SPECIES	Human
ACCESSIONS	BC125039
VECTOR	pDONR223
NCBI SYMBOL	BAG2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcaggcgaagatcaacgctaaagccaacgagggggcgttctccgctcctccatggctgaccgctccagccgctgctggagagcctggaccagctggagctcagggttgaagctttgagagaagcagcaactgctgttgagcaagagagaagaatcctctggaatgat
ccacagtatcctaaatagccagacatgagcagatcagtgacggagagaagaagaatctgactgcaaacgctttgatgggaagaactctcaccgttgaagtgtcagtagaacaattagaaccgccagcagcaagaatccctaaagcatgccacaaggattatgatgaggt
ggcctaataagttctgcatgttgggaatgccaagagtcatttaattgctctacatgcatgttcatctgaggtgccacatggccagttgatcagaagttcaatccatagtaattggctgtgctcttgaagatcagaagaaatgaagagaagattagagactctgcttgaatattgaaa
actctgacaaggccatcaagctattagagcattctaaaggagctgttccaaaactctgcaacaaatgctgaaagcagattcaattag