

## FOXR1-forkhead box R1 Gene

Reference: PTXBC125040

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DLNB13; forkhead box protein R1; forkhead box R1 variant 1; forkhead box R1 variant 2; forkhead box protein N5; forkhead box R1		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	forkhead box R1
<b>GENE ID</b>	283150
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC125040
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	FOXR1
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

atggggaacgagctcttctggccctcaccacatctcactcccccttagcggagcagaaactgccaggtataaactccgaattgtaagccacaaaattacccttagagaaaaacccaaccctgataaggatggtccagattatgagccaacctctggtgtaaatccaacattg  
tgateccccctgaaagctggaggtctcaggacgtaggaagggaggacctgacaagcacactccccctctctcagccacccccagaaggaggaagatgccagctgctcagagccgaggggtggaatcactgccagtcctccagcaagcggctccccctggaagcggttg  
ccctttccccagcactgggagggccccctccagagtcggaggtctcggaagccagcagccagcggggaggctctggtccggccccctcaatctcaactaattgcctggcattaagaaacagttccccctgtggcctcaactgcaacagatctacagttcactcgga  
gctgctggcctgattctgctttgcatgctgtgtccaagccctggccctcagagggagggttcagcttga