



Reference: PTXBC125048

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP10.4; KAP18-4; KRTAP10.4; KRTAP18-4; KRTAP18.4; keratin-associated protein 10-4; high sulfur keratin-associated protein 10.4; keratin-associated protein 18.4; keratin associated protein 10-4		

For research use only.

GENE DESCRIPTION	keratin associated protein 10-4
GENE ID	386672
ORIGIN SPECIES	Human
ACCESSIONS	BC125048
VECTOR	pDONR223
NCBI SYMBOL	KRTAP10-4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr