



Reference: PTXBC125078

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf76; TLT-2; TLT2; dJ238O23.1; trem-like transcript 2 protein; triggering receptor expressed on myeloid cells-like protein 2; triggering receptor expressed on myeloid cells like 2		

For research use only.

GENE DESCRIPTION	triggering receptor expressed on myeloid cells-like 2
GENE ID	79865
ORIGIN SPECIES	Human
ACCESSIONS	BC125078
VECTOR	pDONR223
NCBI SYMBOL	TREML2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]