



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMEM134-transmembrane protein 134 Gene

Reference: PTXBC125134

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

transmembrane protein 134

For research use only.

GENE DESCRIPTION	transmembrane protein 134
GENE ID	80194
ORIGIN SPECIES	Human
ACCESSIONS	BC125134
VECTOR	pDONR223
NCBI SYMBOL	TMEM134
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcgcgccgccgcccacgcttcagcatgatgatgccttcgagctgtccctggaggacggggccctggggcccgagtcacgaggggtgcgcgcgtttggggccgctgcacttcgagcgtcggggccgggttcgaggtggtctgacgaggacaacagcagtcgggtcgtctacacgaa
ctctggagaacgatgatggatggagagccagcctctccgcgagccgatggggagctcgcacacggatcttcacgcgaacttccatccgagctccagctggtctcttcagcaccatcagcagcagcaccagcgtctcttaacaacctgtctgcagctggaccacaacacctttgatcca
gaagaaccgccgagtggtgtcggcctctctctctctctctgtctgggtcgggtgtctcagccagcgtattcttctgtcggggcttctgtgtgtgtcgtcggagtctatcagctgattctctactcgcggctcaaggggccaccggggcttcagttcttactctgccctacttcgaa