

LY6G6D-lymphocyte antigen 6 complex, locus G6D Gene

Reference: PTXBC125138

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf23; G6D; LY6-D; MEGT1; NG25; lymphocyte antigen 6 complex locus protein G6d; lymphocyte antigen-6 G6D; megakaryocyte-enhanced gene transcript 1 protein; protein Ly6-D; lymphocyte antigen 6 complex, locus G6D		

For research use only.

GENE DESCRIPTION	lymphocyte antigen 6 complex, locus G6D
GENE ID	58530
ORIGIN SPECIES	Human
ACCESSIONS	BC125138
VECTOR	pDONR223
NCBI SYMBOL	LY6G6D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaacccagttgttgggatcttctcagctccctgctaggggctgcttggaaacccaatcggtgctacaactgtggtggaagccccagcagttcttgc aaagggccgtgaccacctgtggcgagggcgagacccagccagcctggaacagatcaagctacctggaaccc
ccagtgacctgattcaccaacatccagcctgcctgcagcccatcattgcaatcaagtggagacagagtcggtgggagacgtgactatccagccacagggactgtctacctgggagacctgtgcaacagcgcctgtggcaagccatgtggccctgcaggcatttggctgcagcag
ctaccgccctgacctgtctcttgcaggactgtggagcggatag