

ANP32A-acidic (leucine-rich) nuclear phosphoprotein 32 family, member A Gene

Reference: PTXBC125143

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C15orf1; HPPCn; I1PP2A; LANP; MAPM; PHAP1; PHAPI; PP32; acidic leucine-rich nuclear phosphoprotein 32 family member A; acidic (leucine-rich) nuclear phosphoprotein 32 family, member A; acidic nuclear phosphoprotein pp32; cerebellar leucine rich acidic nuc		

For research use only.

GENE DESCRIPTION	acidic (leucine-rich) nuclear phosphoprotein 32 family, member A
GENE ID	8125
ORIGIN SPECIES	Human
ACCESSIONS	BC125143
VECTOR	pDONR223
NCBI SYMBOL	ANP32A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagatggcgacgagcatttagagctgcggaacagagacccctctgatgtgaaagaactgtcttgacacagtcgggtcgatgaaggcacaactcgaaggcctcacagatgaattgagaactggaattcttaagtacaatcaacgtaggcctcacctcaatcgcaaaactacc
aaagttaacaaacttaagaagcttgaactaagcgataacagagctcagggggcctggaaagtattggcagaaaagtgccgaacctcagcgcataaatttaagtggcaacaaatgaagacctcagcacaatagagccactgaaaaagttagaaaacctcaagagcttagacctttcaa
ttcgaggttaaccaactgaacgactaccgagaaaatgtgtcaagctcctccgcaactcacatattcgcagggctatgacgggacgacaaggaggccctgactcggatctgagggctactgtggaggccctggatgatgaggaggaggatgaggatgaggaggatgatga
agatctcaggtatgtggaagacgaggaggacgaggatgaggaggaggaaagtgaaaggaggacgtgagtgaggaggaggaggatgaagaaggtataacgatggagaggtatgacgaggaagatgaagaagagcttggtgaagaagaagggtcagaagcgaaa
atga