

RP11-297H3.4-hypothetical protein LOC284688 Gene

Reference: PTXBC125144

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC284688		

For research use only.

GENE DESCRIPTION	hypothetical protein LOC284688
GENE ID	284688
ORIGIN SPECIES	Human
ACCESSIONS	BC125144
VECTOR	pDONR223
NCBI SYMBOL	RP11-297H3.4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgattccagatactcaaatggagtgccctctggtacatgagctggatgtgctggcagcttccatctggctggcagctatcctgaggttagcaaaagcagaaggcaaatcttggggagctaaaagacatgtgcacgtgccagttttgacaaatgctcagtttaattgtttctatgattctt
ctcccagagacacaaactgtgtgcaccaggcaccagaactatccccatccccaccacccagatcccaagggaacttcataattacacaggtttggccaacatctaccccacaggaaatggacacatcagcccctaggatgaagaatgatgcagggcagggaaagcagagaaa
tctcctgagtgagaagaagacaaattcacagaagcctga