

STELLAR-germ and embryonic stem cell enriched protein STELLA Gene

Reference: PTXBC125145

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	germ and embryonic stem cell enriched protein STELLA
GENE ID	400206
ORIGIN SPECIES	Human
ACCESSIONS	BC125145
VECTOR	pDONR223
NCBI SYMBOL	STELLAR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaccatcacagttaatccaactacatccagggtctccacaaatgtcacgaagaaatccgggacgattcaggggcctctcaaatctctccgagacgttgataaagaaccttagtaacttgactatcaacgctagtagcgaatctgtttccctctattggaagctttactccgtc
gagagctctgtggggcagcagctcctcagggaatcgaagatgagtgctttacagcaggagaggagtaagaacactgctgtctgtgcagagagaaaaagatggcaagattgagatacatgttactggcgaggttcgtacgcatgaaagaagaccaacaacaaggagcctaaggag
ttaagaaggaaatcaagaccattcaaatgtccctgcagttctgcgtgtctaatggatgggaccttctgagaatgctagaatagagaatcaagacaccaagccactcagccataa