



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SLC38A11-solute carrier family 38, member 11 Gene

Reference: PTXBC125157

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AVT2; solute carrier family 38 member 11		
For research use only.			

GENE DESCRIPTION	solute carrier family 38, member 11
GENE ID	151258
ORIGIN SPECIES	Human
ACCESSIONS	BC125157
VECTOR	pDONR223
NCBI SYMBOL	SLC38A11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr