

KRT8P12-keratin 8 pseudogene 12 Gene

Reference: PTXBC125159

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	keratin 8 pseudogene 12
GENE ID	90133
ORIGIN SPECIES	Human
ACCESSIONS	BC125159
VECTOR	pDONR223
NCBI SYMBOL	KRT8P12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgaggcatcaccacatcacggttaaccagagcctgtgagccccccttaacctggaggtggactccaacattcaggctgtataccaggagagaggagcagatcaagaccctcaacaacaagtttgccctctcacagaccaggtatggttcctggagcagcagaacaaggtgctggagaccaagtggagcctctgcagcagcacaagatggctgggagcaacatggacaacatgttcggagctatatcaacaaccttaggcggctgctggagccaggagaagctgaggctggaagcggagcttggcaacatgcaggggctggtggaggtcagggaacaagtatga