

## IL22RA2-interleukin 22 receptor, alpha 2 Gene

Reference: PTXBC125167

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRF2-10; CRF2-S1; CRF2X; IL-22BP; IL-22R-alpha-2; IL-22RA2; ZCYTOR16; interleukin-22 receptor subunit alpha-2; cytokine receptor class-II member 10; cytokine receptor family type 2, soluble 1; interleukin 22 receptor, alpha 2; interleukin 22-binding prote		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	interleukin 22 receptor, alpha 2
<b>GENE ID</b>	116379
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC125167
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	IL22RA2
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon )

atgatcctaaacattgcttctagcttctcctcagtttctcttactggtgtagcaggaactcagtcaacgcagtgctctggaagcctcagagggtacaattcagtcgccgaaatttcacaacatttgcattggcagcctgggagggcacttactggcaacagcagtgctctatttgtgcag  
tacaaaaatatafgacagagacaatggaaaaataaagaagactgttggggtactcaagaactctctgtgaccttaccagtgaaacctcagacatcacaggaaaccttattacggagggtgagggcgccctcgctgggagctactcagaatggagcatgacgccggttccctcctgt  
gggaacaaaaatagatcctccagtcataataaccaagtcattggctcttctgtgtaattctccatgctcctcaatttaccatagatacacaaggaaaaaatgtatctatagaagattactatgaaactactataccagagttttataaatacaattcactagaaaaggagcaaaaggttta  
tgaagggtctcacagagcgttgaattgaagcttaaccacactccagctactgtgtagtgctgaaatatacagccatgttagacagaagaagtcagagaagtgaagagagatgtgtgaaattccatga