

C13orf37-chromosome 13 open reading frame 37 Gene

Reference: PTXBC125183

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C13orf37; MOZART1; mitotic-spindle organizing protein 1; mitotic-spindle organizing protein associated with a ring of gamma-tubulin 1; mitotic spindle organizing protein 1		

For research use only.

GENE DESCRIPTION	chromosome 13 open reading frame 37
GENE ID	440145
ORIGIN SPECIES	Human
ACCESSIONS	BC125183
VECTOR	pDONR223
NCBI SYMBOL	C13orf37
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgagtagcggcggtgctggggcggcgccggcgccggcgccggaatctgaatgcggcgccggcgagaccatggacgttctgctgagattcaagaatttgaatactggcttagatatgaaactctgtctattgtgtacggctttgtgaacaagggaattaccagaagctttatca
tcggttattaaaggagcttcgaaggctactgaagcactgaaggctgctgaaaatatgacaagctga