



Reference: PTXBC125184

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IRAK-2; interleukin-1 receptor-associated kinase-like 2; interleukin 1 receptor associated kinase 2		

For research use only.

GENE DESCRIPTION	interleukin-1 receptor-associated kinase 2
GENE ID	3656
ORIGIN SPECIES	Human
ACCESSIONS	BC125184
VECTOR	pDONR223
NCBI SYMBOL	IRAK2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]