

RACGAP1P-Rac GTPase activating protein 1 pseudogene Gene

Reference: PTXBC125190

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	Rac GTPase activating protein 1 pseudogene
GENE ID	83956
ORIGIN SPECIES	Human
ACCESSIONS	BC125190
VECTOR	pDONR223
NCBI SYMBOL	RACGAP1P
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagcagcagaatcacagatgaagacaacagcatatctgcatgtaccaggtgtgtgtaactgccccagggccaacagggaacattggtttctcatgattcacttgcagagagtggtcagagtcataactataaaatgaatgttgccaatctggctgaagtctttggcttacaat
agtggcccatgtgtgccaatccagaaccagtgacgatgttacaggacatcaagtgtcaaccaaggtgtgagcgcctgcttctgccccctggagtactggagtcagttctcaatagtgagcagagaaacattgacccccctacatgtcactgaaaactcaaatgccttttcaacacca
gatattaaagtga