



Reference: PTXBC125274

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ANKRD41; LEM3; LEMD6; ankyrin repeat and LEM domain-containing protein 1; LEM domain containing 6; LEM-domain containing 3; ankyrin repeat domain-containing protein 41; ankyrin repeat and LEM domain containing 1		

For research use only.

GENE DESCRIPTION	ankyrin repeat and LEM domain containing 1
GENE ID	126549
ORIGIN SPECIES	Human
ACCESSIONS	BC125274
VECTOR	pDONR223
NCBI SYMBOL	ANKLE1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]