



Reference: PTXBC126130

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COE3; EBF-3; OE-2; transcription factor COE3; olf-1/EBF-like 2; early B-cell factor 3		
For research use only.			

GENE DESCRIPTION	early B-cell factor 3
GENE ID	253738
ORIGIN SPECIES	Human
ACCESSIONS	BC126130
VECTOR	pDONR223
NCBI SYMBOL	EBF3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr