



Reference: PTXBC126141

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AFF29; CDC14B; zinc finger protein 367; C2H2 zinc finger protein ZFF29; C2H2 zinc finger protein ZFF29a; CDC14 cell division cycle 14 homolog B		

For research use only.

GENE DESCRIPTION	zinc finger protein 367
GENE ID	195828
ORIGIN SPECIES	Human
ACCESSIONS	BC126141
VECTOR	pDONR223
NCBI SYMBOL	ZNF367
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]