

C7orf53-chromosome 7 open reading frame 53 Gene

Reference: PTXBC126149

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C7orf53; leucine-rich single-pass membrane protein 1; leucine rich single-pass membrane protein 1		

For research use only.

GENE DESCRIPTION	chromosome 7 open reading frame 53
GENE ID	286006
ORIGIN SPECIES	Human
ACCESSIONS	BC126149
VECTOR	pDONR223
NCBI SYMBOL	C7orf53
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactcattctccagacactggttctgtggcattcaggaagatggaagctttatgtgtggtgattccataatgacttaacaaactaaactctgtccagccgcatcgagcatctgttcctctagaggacaaaatccagtccttggcacaactcaggaaatggaagccggagtctgtttttgtgggctgctaattgtgctgattgtcagcctggcactggttttttcgtgataatttctaatagttcaaaactggaacaagatggatgatgtgtcaagaagactaacagctgaaggaaaagacatagatgatcttaagagaatcaataacatgatcgtaagcgactcaaccaactcaacaactggactctgaacaaaactaa