



Reference: PTXBC126170

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	prolyl 4-hydroxylase subunit alpha-3; 4-PH alpha-3; C-P4H alpha III; collagen prolyl 4-hydroxylase alpha(III); procollagen-proline, 2-oxoglutarate 4-dioxygenase (proline 4-hydroxylase), alpha polypeptide III; prolyl 4-hydroxylase, alpha polypeptide III; p		

For research use only.

GENE DESCRIPTION	prolyl 4-hydroxylase, alpha polypeptide III
GENE ID	283208
ORIGIN SPECIES	Human
ACCESSIONS	BC126170
VECTOR	pDONR223
NCBI SYMBOL	P4HA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]