

C21orf99-chromosome 21 open reading frame 99 Gene

Reference: PTXBC126176

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	chromosome 21 open reading frame 99
GENE ID	149992
ORIGIN SPECIES	Human
ACCESSIONS	BC126176
VECTOR	pDONR223
NCBI SYMBOL	C21orf99
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaaatgcaaatggagttgataagtttaaatgcatcatcaacaacttttggatataaacaagaatatactaaaaattctcaaatagtaatccagaaggacacatctgaaggacacctgacgaggctgcacccctggcggaaagaacacctgacacggctgaaagcttggtggaagaa
cacctgacgaataggatagcagtgga