

## IL29-interleukin 29 (interferon, lambda 1) Gene

Reference: PTXBC126183

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL29; IL-29; interferon lambda-1; IFN-lambda-1; cytokine Zcyto21; interleukin 29 (interferon, lambda 1); interleukin-29; interferon lambda 1		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	interleukin 29 (interferon, lambda 1)
<b>GENE ID</b>	282618
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC126183
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	IL29
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

atggctcagctggaccgtggtgctggtgactftggtgctagcttggccgtggcagccctgtcccactccaagcccacacaactgggaaggctgccacattggcaggttcaatctctgtcaccacaggagctagcagcttcaagaaggccaggagccttgaagagtca  
ctcaagctgaaaaactggagttgacgctctctgtctccccgggaattgggacctgaggtctctccaggtgaggagcgcctgtggccttggaggctgagctggccctgacgctgaaagtcctggaggccgctgctggccagcctggaggactcttagaccagccctcacac  
cctgaccacatctctccagctccaggcctgtatccagctcagcccacagcagggccagcccccggggccgctccaccactggctgcaccggctccaggagcccccaaaaggagtcgctgctggcctggaggatctgtcacctcaacctctccctctcacgcga  
gacctcaaatatgtgccgatgggacctgtgtctgagaagctcaaccacctgagtcacctga