

C3orf55-chromosome 3 open reading frame 55 Gene

Reference: PTXBC126187

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C3orf55; PQ-loop repeat-containing protein 2-like; PQ loop repeat containing 2 like		
For research use only.			

GENE DESCRIPTION	chromosome 3 open reading frame 55
GENE ID	152078
ORIGIN SPECIES	Human
ACCESSIONS	BC126187
VECTOR	pDONR223
NCBI SYMBOL	C3orf55
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagtcgtgggaaactatagagtaacactgcaaaactcaagcactgatactgccggagagcactgacctgccttagaagtcagctctttagcctacagaaatggaagagtgatgaagcagtcctctgggtttctggattgctggatagggtgagacctgacaaattcaagctgctacctgactaaccaactgcctattcagattttacagccatcttcgacatgaacacggatgtaatcactctcacaattcatgtactacaggttaaagaatcagaagaaaaattctgaagaagatctacggaaggaactcatttctgctgtttgcatggccatga