

C21orf94-chromosome 21 open reading frame 94 Gene

Reference: PTXBC126226

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf94; NCRNA00314; long intergenic non-protein coding RNA 314		
For research use only.			

GENE DESCRIPTION	chromosome 21 open reading frame 94
GENE ID	246705
ORIGIN SPECIES	Human
ACCESSIONS	BC126226
VECTOR	pDONR223
NCBI SYMBOL	C21orf94
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgaagattgcagtgatggaagtgtctcagtcagcctggaagactgtatggaaggaggtgcctggccacactgtcatgtccctctccctaaaaaggacacatttgcccttgcaggcaacacatacatttggtctgtgtatatacttggctctgtattctgtatttaacattctaa