

LOC339524-hypothetical protein LOC339524 Gene

Reference: PTXBC126240

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC339524		
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC339524
GENE ID	339524
ORIGIN SPECIES	Human
ACCESSIONS	BC126240
VECTOR	pDONR223
NCBI SYMBOL	LOC339524
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgccccgccgacgtgggaactgtgccacatggagtgctgggagctcattgccgcgagttctacaccacagggccaggctgttgcctcccatcggtcgtgccccagcacctgtgttattaaggactcattgtgagcgccatcattacaagggtgtgggtactacat
atactccctattttctatcttgcgaaggctgcagggcgcatgcacagtgctgtgcacgggtgggctagtgctagcccgaggagcgagcgggggccggcagggcggtgggcgcggcctcgtctcgtgctgctgagcgtgtccggcgccggcgaggccct
ggcgtcttcacccctggagggcggtggggcggtgctggggcgcttccctcctcagggccctcaaacgcgaaggggtttccgcttccagtcgatggtcgtctcctccatccccggctgcatctccattaccgctgccatttcttggccatccaggctcttggctccactg
gggtctcgttcccttcccggttccccctccaggtcggtccttcttgcaggactacgcaggggcttgacccagggcgctgttttaggcggatctgggtcccttgcactccagggttcttccacttccgaattctggagaaccgggaatcaagccctgcggtctcttcttct
ccttctgcccgaagcacgctcatgtctgcagggcacagttctgaaagtgcgggagagaacaaggagttcttcttctccagagttgcttctttaa