

CCL8-chemokine (C-C motif) ligand 8 Gene

Reference: PTXBC126242

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HC14; MCP-2; MCP2; SCYA10; SCYA8; C-C motif chemokine 8; chemokine (C-C motif) ligand 8; monocyte chemoattractant protein 2; monocyte chemotactic protein 2; small inducible cytokine subfamily A (Cys-Cys), member 8 (monocyte chemotactic protein 2); small-i		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 8
GENE ID	6355
ORIGIN SPECIES	Human
ACCESSIONS	BC126242
VECTOR	pDONR223
NCBI SYMBOL	CCL8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggtttctgcagcgttctgtgctgtctcatggcagccatttcagccctcaggagctgtctcagccagattcagttccattccaatcacctgtctttaacgtgatcaataggaaaattctatccagaggctggagagctacacaagaatcaccaacatccaatgtcccaaggaa
gctgtgatcttcaagccaacggggcaaggaggtctgtgtgacccaaggagagatgggtcagggaattccatgaagcatctggaccaaatatttcaaaatctgaagccatga