

TMIE-transmembrane inner ear Gene

Reference: PTXBC126258

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DFNB6; transmembrane inner ear expressed protein; transmembrane inner ear protein; transmembrane inner ear		

For research use only.

GENE DESCRIPTION	transmembrane inner ear
GENE ID	259236
ORIGIN SPECIES	Human
ACCESSIONS	BC126258
VECTOR	pDONR223
NCBI SYMBOL	TMIE
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggggtggcgggcgccgggtccctctgctgctgggcgccgacctgggggtgctcctcgggggtgccgggcagctggtggagcccagcacggccccccaagcccaagccctcctcctgaccaaggagacagtgggtgttctgggacatgcgctgtggcac
gtggtggcatctttcgtcttcgtgttgcacatcatcacgctgtgctgtgtctcaactgctgtgcccaggaccggaaaggagatcgaagcccgtacctgcagcgaagcagccaagatgtacacagacaagctggagactgtcccaccctcaatgagctcacagaagtccc
aggagaggataagaagaagaagaagaagaaggacagctgtggacacagtgccatcaaagttagaggagatgagaagaatgaggccaagaagaagaaggagagaatga