



Reference: PTXBC126259

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EEF-2; EF-2; EF2; SCA26; elongation factor 2; polypeptidyl-tRNA translocase; eukaryotic translation elongation factor 2		

For research use only.

GENE DESCRIPTION	eukaryotic translation elongation factor 2
GENE ID	1938
ORIGIN SPECIES	Human
ACCESSIONS	BC126259
VECTOR	pDONR223
NCBI SYMBOL	EEF2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]