

C22orf31-chromosome 22 open reading frame 31 Gene

Reference: PTXBC126264

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C22orf31; HS747E2A; bK747E2.1; chromosome 22 open reading frame 31		

For research use only.

GENE DESCRIPTION	chromosome 22 open reading frame 31
GENE ID	25770
ORIGIN SPECIES	Human
ACCESSIONS	BC126264
VECTOR	pDONR223
NCBI SYMBOL	C22orf31
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacccaatcaatgtgagacgagacccagcatccctctatgtgactccgacagtcctcttatttaataaccaggctcagagactgctatgtggaactaccggctctcaccaacatctggatggccagaacatgtgcaagcagaacattatgccccagcaccagctaccactctctt
gggaagtgttaaggaaacccatttaattgccagttctctcctgtgtaagctgtactcagcgccgaactgaagaataagtgctgccaccaccatgcaagtttgagaaggaaaactctcgaagagattaaagcacaaaggacgattcagtgatgaaagccaccagcaggccaggaaaag
aaacttcacagttccaagagcaagcagccagcagggcatagggagcctgcaggagcctcagagagagtaaaagaaagtcaaggagaaaaaactaacagtcgcgaagatcttgaggacagatatgctgaacatgtggctgccaccaagcgctacccagcagtgaggacag
cagcctggaaggccgagtggtgcttctgaaacccaaagagacagcagttgtcggaggacacgctaaccatcctggtctccccacagagggttaccaggctctgtaccacgctgtgtggagccaatgctgtggaatcctcagggaccccccaagaggtacagcctggagctggg
caaggccattaaacaaaagctctggaggctcttgcagtcaggggtccatctgaagggtgctcagagggaccggttccctggcaggaagcagccagggtgtccacagaggagcctgtactcaagaatggccaagttaaagagcaaaaatag