

HSPB3-heat shock 27kDa protein 3 Gene

Reference: PTXBC126272

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DHMN2C; HMN2C; HSPL27; heat shock protein beta-3; HSP 17; heat shock 17 kDa protein; heat shock 27kDa protein 3; protein 3; heat shock protein family B (small) member 3		

For research use only.

GENE DESCRIPTION	heat shock 27kDa protein 3
GENE ID	8988
ORIGIN SPECIES	Human
ACCESSIONS	BC126272
VECTOR	pDONR223
NCBI SYMBOL	HSPB3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaaaatcatttggagcaccatagagattccagtgcttaccaggaaggttgaagctcgaggtctagaagactgcaggtggtatcatgctttatatgcactgctgggccaaccatcgtggacctgaggaaaaccaggcgagcgcagtcctccagtggtactcagcggcag
agacgccaccgccgagaaggcaaatccactttcagatcctgctggacgtggtccagttcctccctgaagacatcatcattcagaccttcgaaggctggctactgataaaagcacaacacggaaccagaatggatgagcacgggtttatctcaagaagcttcacccgacagtacaaactacca
gatggcgtggaaatcaagatttgcagtcctcgtcatgatggaatttgggtgggaagtaaggatccagttgggactaagtga