

CLEC12A-C-type lectin domain family 12, member A Gene

Reference: PTXBC126289

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD371; CLL-1; CLL1; DCAL-2; M1CL; C-type lectin domain family 12 member A; C-type lectin protein CLL-1; C-type lectin-like molecule-1; dendritic cell-associated lectin 2; myeloid inhibitory C-type lectin-like receptor		

For research use only.

GENE DESCRIPTION	C-type lectin domain family 12, member A
GENE ID	160364
ORIGIN SPECIES	Human
ACCESSIONS	BC126289
VECTOR	pDONR223
NCBI SYMBOL	CLEC12A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgaagaagtacttatgcagatctcaattccagaactccagtgagatggaaaaatccagaaattggcaaatggggaaaaagcactccagctccctctcatgtatggcgtccagcagcctgtttctgactcttctgtgcctctgtgctcattggattggagcttggcaagcat
gtttcacgtaactttgaagatagaaatgaaaaaatgaacaaactacaaaacatcagtgagagctccagagaaatatttctacaaactgatgagtaacatgaatctccaacaagatcaggaaacctctccaccacactgaaacaatagccacaaattatgtcgtgagctatagcaaag
aacaagagcacaatgtaagcctgtccaaggagatggaattggcataaggacagctgtatttcttaagtgtatgtccaacatggcaggagagtaaaatggcctgtgtcgtcagaatgccagcctgttgaagataaacaacaaaatgcaattggaattataaaatccagagtagatc
atafgactattggctgggattatctcctgaagaagattccactcgtggatgagagtgataataatacaactcctctgcctgggtataagaacgcacactgacttaataacatgtattgtgataataaatagactatattgtcaattatcactgcactataaacaagaatgatattgaga
agatggccaatccagtcagctgtgtctacatatttagggagcatga