

PI15-peptidase inhibitor 15 Gene

Reference: PTXBC126290

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRISP8; P24TI; P25TI; peptidase inhibitor 15; 25 kDa trypsin inhibitor; CRISP-8; PI-15; cysteine-rich secretory protein 8; protease inhibitor 15; sugarCrisp		

For research use only.

GENE DESCRIPTION	peptidase inhibitor 15
GENE ID	51050
ORIGIN SPECIES	Human
ACCESSIONS	BC126290
VECTOR	pDONR223
NCBI SYMBOL	PI15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgataagcaatctcgcgtcagcagtgcaactctgttctccctctctgtgaagcaagtacgtctcctcaattccactgactcatccccgccaaccaataattcactgatattgaagcagctctgaaagcacaattagattcagcggatccccaaagccaggcggaagcgctacatt
cgagaatgacatgatcgccattcttgattatcataatcaagttcgggcaaaagtgttccaccggcagcaaatatggaatatatgtttgggatgaaatcttgcacaaatcggcagagccttggcgctacttgcaattggaccatggaccttctacttactgagatttttggccaaatct
atctgtacgactggaagatactctctatttccagttggtaagccatggatgatgaagtgaagattatgctttccatccccaggattgcaacccagatgtcctatgagatgttttggccatgtgcacacattatagcagatgtttggccacttccaatcggataggatgcgcaa
ttcatacttgccaaaacatgaatgtttgggagctgtgtggcgacgtgcagtttacttggtatgcaactatgcccaaaaggcgaattggattggagaagcaccatataaagtaggggtaccatgttcactgtctccaagttaggggatctgtactgacaatctgtgtttccaggagttacgt
caaacactctgtactggtttaataaa