



Reference: PTXBC126307

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MYCD		
For research use only.			

GENE DESCRIPTION	myocardin
GENE ID	93649
ORIGIN SPECIES	Human
ACCESSIONS	BC126307
VECTOR	pDONR223
NCBI SYMBOL	MYOCD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]